

Advances In Prosthodontics™

New Findings & Best Practices for Cosmetic & Restorative Dentistry

INSIDE THIS ISSUE

CASE SUMMARY:

Implant Overdentures after Failed Implant-Fixed Prosthesis

Background
pages 1-2

Treatment Benefits,
and Final Result
pages 2

PRODUCT SPOTLIGHT:

MAGNET-X™
Implant Overdenture
page 3

EXCITING NEWS:

Coming this Fall...
New Office Location +
Updated Name
page 4

© Sams Media Group, LLC

Do You Have a Patient with Failing Restorations?

MAGNET-X
Implant
Overdenture



Final Result

Here's the story of a patient who was referred from another dental professional...

If you have a patient with several missing, broken or severely worn teeth — or failed restorations — they may be a candidate for oral rehabilitation with the assistance of an advanced Prosthodontist.

Read more on the next page...

◆ Do you want another opinion on a complex case? Or to talk and share ideas?

There are many ways we can help each other by working together on complex cases and by sharing knowledge and experience.

If you'd like to discuss a treatment, or talk on the phone about how we can help each other, please call or email our office with specific days and times when you're available.

More Examples

Stay informed on new findings and best practices for cosmetic and restorative dentistry.

See more examples of complex cases online:



AdvancedProsthodontist.com
/ PleasantonPros

Implant Overdentures after Failed Implant-Fixed Prosthesis

by the Specialists of *Pleasanton Prosthodontics* — Dr. Maggie Chao & Dr. Olivia Nguyen

Patient “Sharon” was 62-year-old when she came to my office with a failed implant-fixed prosthesis. Due to a history of peri-implantitis, we elected to restore her with a **MAGNET-X™ Implant Overdenture**.

Implant Overdentures are an excellent solution for patients with prior failed implants due to their superior cleanability and efficient force distribution.

- **Cleanability due to removability** — Unlike fixed dentures, overdentures can be easily removed for thorough cleaning, ensuring better oral hygiene and reducing the risk of peri-implant diseases.
- **Balanced force distribution** — implant overdentures distribute occlusal forces more evenly across the jawbone than a fixed prosthesis, which (1) helps maintain the structural integrity of the implants and surrounding bone, and (2) minimizes risk of overloading by any single implant, thus increasing the overall longevity.

If you have patients wanting a reliable and maintainable restorative solution, I suggest looking into Implant Overdentures. We can work with you or assist your patients with this treatment.



Final Result



We give all referred patients a warranty and refer them back to the general dentist for regular hygiene unless advised otherwise by the referring office, depending on the patient and treatment. If you have questions about our prosthodontic protocols, please email us at: office@PleasantonPros.com

◆ Do you have a patient in Sharon’s situation? Or another complex case?

If you have a question about a case and want another opinion, give our office a call. Or, if there’s a case that you want to work together on, please fill out and send in the enclosed referral form.

We will take great care of your patient and keep you informed. Dental professionals refer to us as **The Complex Case Specialists™** because we perform complex cases every day.

MAGNET-X™ Implant Overdenture

Source: pantheradental.com



Using conical magnets in its patent-pending system, the Magnet-X™ offers stable forces, lateral guidance and stability, and is self-aligning to make it easy to handle for patients with less dexterity. Here's how it works...


Magnetic

- The magnetic forces guide the denture to the right position and keep it in place. The magnets are self-aligning the bar into right position autonomously. The patient does not need to find the right position or to chew the denture in place. The conical magnetic male and female parts provide a durable solution, constant retention and lateral guidance and stability.



Custom

- Directly on the implants, the custom abutments adjust the angle of the implants up to 30° in order to make the whole restoration perpendicular to the occlusal plane. This is the role of this first piece, which can be adapted to a wide range of implant connections. Since the custom abutments are straight on the implants, you do not need additional Multi-Units.



Durable

- The Magnet-X™ implant bar includes durable titanium magnet caps and attachments from German company steco-system-technik for a constant retention with no wear parts to replace. Due to the titanium housing, the magnets are corrosion proof. This is a world's first that respond to geriatric patient needs. Every detail has been taken into consideration to create an impeccable product.



With 2 prosthodontists in the same specialty practice, we have increased capacity to work with you on challenging treatments, so your patients don't have to wait weeks to get a plan for comprehensive care and/or to address reconstructive emergencies. In addition, we have digital workflows that reduce the number of appointments even for comprehensive treatment, while still maximizing function and esthetic outcomes, so your patients can receive restorative care in fewer appointments.

If you have questions about our prosthodontic protocols, please email us: office@pleasantonpros.com



◆ Why other dental professionals work with a Certified Prosthodontist

Certified Prosthodontists are specialists in implant, esthetic, and reconstructive dentistry.

1. You have a patient requiring treatment outside your typical scope of practice.
2. You have a patient with complex needs that may drain your time.
3. You have a demanding patient who wants perfect esthetics.
4. You have a question and want to discuss a case with a colleague to ease your mind.

New Office Location + Updated Name

Same great specialty service, always served with a smile...
We are looking forward to better helping you and
your referred patients at our expanded facility!



Dr. Maggie Chao and Dr. Olivia Nguyen at **Pleasanton Prosthodontics (soon to be Renu Prosthodontics)** are specialists in prosthodontics and are available to help you and your patients with complex cases (yes, even during this office relocation!).



PLEASANTON
PROSTHODONTICS
CENTER FOR ADVANCED DENTISTRY

1475 Cedarwood Lane, Suite C
Pleasanton, CA 94566

Coming
this
Fall

RENU
PROSTHODONTICS

2600 Old Crow Canyon Rd # 100
San Ramon, CA 94583

Coming Soon — full announcement with photos of our new office building!

Please call our office at (925) 425-7545 to refer a patient or with any questions.

◆ How can a Prosthodontist assist you and your dental team?

We specialize in the treatment of complex cases, usually involving several procedures over months of care. The next time you see a challenging case, please feel free to call us and we can discuss treatment planning or help you with any part of the treatment. Our goal is to be a resource for your office by helping manage difficult and time-consuming procedures, restoring the function and esthetics that the patient desires, then referring them back to your office for their ongoing care.